



## Smoking kills

### Expired mode:



Awareness for smoking cessation.  
Monitoring during withdrawal programmes.  
Helps motivation.  
Awareness for passive smoking.

Screening for accidental intoxication.

### Monitoring mode:



Awareness in public places for the respect of smoking/non-smoking areas.

In monitoring mode, the Tabataba™ records measured CO values within an adjustable time range, thus giving a good indication of induced passive smoking.

### Alarm mode:



In alarm mode, the Tabataba™ sets off a visual and sound alarm to warn the wearer of the device that the preset threshold has been reached.



## CO Tester for smoking awareness

Real size

# Tabataba™

## Features



- 3 modes: Expired, Monitoring, Alarm
- CO Measurement: 0 - 400 ppm
- Expired Mode Calculation: %CoHb, %CoHbfetal
- Monitoring Mode Calculation: TWA, STEL, C
- Resolution: 1 ppm
- Weight: 6 oz
- Size: 4.33" x 2.75" x 1.18"
- Cell type: Electrochemical
- Usual life span of cell: 3 years (depending on use)
- Measurement time: 15 seconds
- Screen: LCD
- Batteries: 2 x AA 1,5 volts LR6 Alkaline
- Communication: Serial port (software and lead not included)
- Parametering: Date, time
- Recording times: 1/4hr, 1hr, 8hr, 24hr
- Adjustable alarm thresholds

DEPISTEO LCC  
3495 Piedmont Rd. NE  
Bldg 11, Suite 710  
Atlanta, GA 30305 - USA  
Tel. (404) 495-5911 - Fax (404) 233-8625  
[www.despisteo.com](http://www.despisteo.com)



DEPCTAB\_001